

Your
Vision
in to the
Future



Light Weight



Flexible & Compact



User-friendly interface



Simplify procedures with
smooth, safe integration.



Certified & Approved By:




Technical Specifications

Specifications	High Frequency Mobile X-Ray Generator		
Model	EP Corsa 4 SM	EP Corsa 6 SR	EP Corsa 6 SM
Application	General Radiography		
Max Output (kW)	4	6	6
kV Range.	40 - 120	40 - 120	40 - 110
mA Max.	100	110	150
Digital Display	kV, mA & mAs		
X-Ray Tube	Stationary Anode		
Control panel	Feather Touch (Micro Processor Controller Based)		
Collimator	Manually / Motorized Light Beam Diaphragm (LBD)		
APR Facility	Available		
Power Required	1Ø 240 VAC±10% 50 Hz with 15 Amps Power Socket (Required Proper Line Earthing)		
Tube Stands	Gas Spring Balancer Based		

Due to continuous R & D or product development and improvement design, specifications, appearance & Data indicated may vary from the product illustrated in the catalogue.

Clinical Images



Epsilon Healthcare Solutions (P) Ltd.

Corporate Office : E.D Tower, F-101, J-17 D, Ramesh Marg, C-Scheme, Jaipur-302001, Raj. (INDIA) Tel : +91-141-4910019

Registered Office : G-1-57, RIICO Industrial Area Extension, Phase-2, Bagru, Jaipur, Raj. (INDIA)

Works : G-1-58, RIICO Industrial Area Extension, Phase-2, Bagru, Jaipur, Raj. (INDIA)

E-mail : marketing@epsilonhealthcare.com, **Web** : www.epsilonhealthcare.com

Marketed by:

Toll Free : 1800 1234 122