

passion to lead
EPSILON



EP Corsa Economical Mobile DR System

- EP Corsa Mobile digital Radiography system is a high-frequency Radiography solution ideal for emergency rooms, intensive Care Units (ICUs), and bedside X-Ray procedures.
- It delivers high-quality diagnostic images while ensuring minimal radiation exposure, prioritizing patient safety and comfort.
- The system's mobility allows healthcare professionals to efficiently perform xrays in various hospital environments, enhancing patient care in critical and urgent situations.

Why Epsilon !



Clinical Images

Our X-ray machines are engineered to deliver high-resolution, detailed images that reveal even the most minute anatomical structures. This level of clarity is crucial for accurate diagnosis and effective treatment planning, ultimately improving patient care.



Hand AP



Cervical Spine



Chest PA



LS Spine Lateral



Knee AP




KUB

Application Feild



Smooth execution



Technical Specifications

Specifications	Mobile X-Ray Generator		
Model	EP Corsa 4 SM	EP Corsa 6 SR	EP Corsa 6 SM
Max Output (kW)	4	6	6
kV Range	40-120	40-120	40-110
mA Max.	100	110	150
X-Ray Tube Head	Monoblock Version		
X-Ray Tube	Stationary Anode		
Tube Stand	Gas Spring Balancer Based		
Digital Display	KV, mA & mAs		
Control Panel	Feather Touch (Micro Processor Controller Based)		
Imaging Console	Touch Screen		
APR Facility	Available		
Collimator	Manually / Motorized Light Beam Diaphragm (LBD)		
Power Required	1Ø 240 VAC±10% 50 Hz with 15 Amps Power Socket (Required Proper Line Earthing)		
	Flat Panel Detector		
Scintillator	CsI (Cesium Iodide)		
Detector Size	14x17 Inch / 17x17 Inch		
Type	Wireless		

Overview



Due to continuous R & D or product development and improvement design, specifications, appearance & Data indicated may vary from the product illustrated in the catalogue.

Epsilon Healthcare Solutions (P) Ltd.

Corporate Office : E.D Tower, F-101, J-17 D, Ramesh Marg, C-Scheme, Jaipur-302001, Raj. (INDIA) Tel : +91-141-4910019

Registered Office : G-1-57, RIICO Industrial Area Extension, Phase-2, Bagru, Jaipur, Raj. (INDIA)

Works : G-1-58, RIICO Industrial Area Extension, Phase-2, Bagru, Jaipur, Raj. (INDIA)
E-mail : marketing@epsilonhealthcare.com, Web : www.epsilonhealthcare.com

Auth. Channel Partner

Toll Free : 1800 1234 122