



INTRODUCTION

PINAKI SURGICAL IS ONE OF THE PIONEERING MEDICAL EQUIPMENT COMPANIES IN INDIA, TO PROVIDE HIGH QUALITY, AFFORDABLE NEONATOLOGY EQUIPMENTS.

SINCE 2014, WE HAVE BUILT UP AN IMPRESSIVE ARRAY OF PRODUCTS FOCUSING PRIMARILY ON MEETING THE NEEDS OF NICU'S ACROSS THE COUNTRY. OUR SUCCESS HAS BEEN HARD EARNED. WE HAVE WORKED HARD TO CONTINUOUSLY IMPROVE OUR PRODUCTS, BY LISTENING TO OUR CUSTOMERS, UNDERSTANDING THEIR UNIQUE REQUIREMENTS AND TRYING TO REPLICATE THOSE NEEDS.

OUR LED BASED PHOTOTHERAPIES ARE RATED AMONG THE BEST PERFORMING MACHINES IN INDIA, PROVIDING RELIEF TO BABIES WITH SWIFT & EFFECTIVE PERFORMANCE. OUR PRODUCT FACILITIES MEET INTERNATIONAL STANDARDS AND OUR PRODUCTS ARE CERTIFIED BY NATIONAL & INTERNATIONAL CERTIFICATION AGENCIES.

PINAKI SURGICAL IS ALSO ONE OF THE LEADING MANUFACTURERS OF INFANT RADIANT WARMERS WITH INSTALLATIONS ACROSS REPUTED HOSPITALS ALL OVER INDIA. IN ADDITION TO OUR EQUIPMENT RANGE, WE STARTED TO PROVIDING SIMPLE YET VERY EFFECTIVE SOLUTIONS FOR NEONATOLOGISTS & NURSING STAFF BY LAUNCHING OUR NICU VITALS RANGE.



BABY WARMER

MODEL: ELITE



- Ceramic Heater
- 183cm (H) x 65cm (W) x 88cm (L)
- Easy mobility over a bed. Bassinet or examination table.
- Microprocessor Controlled
- Set Temperature range 32°C to 39.9°C.
- Temperature Alarms $\pm 0.5^{\circ}\text{C}$.
- Alarms- High/Low Temperature, Sensor Fail, Power Fail (Servo Mode)
- Programmable Features – Set Temperature (High/Low Limit),

BABY WARMER

MODEL: ECONOMY



- Ceramic Heater
- 183cm (H) x 77cm (W) x 100cm (L)
- Easy mobility over a bed. Bassinet or examination table.
- Microprocessor Controlled
- Set Temperature range 32°C to 39.9°C.
- Temperature Alarms $\pm 0.5^{\circ}\text{C}$.
- Alarms- High/Low Temperature, Sensor Fail, Power Fail (Servo Mode)
- Programmable Features – Set Temperature (High/Low Limit), Mute, Manual Time, Apgar Time.

BABY WARMER

MODEL: DELUXE



- Ceramic Heater
- 183cm (H) x 77cm (W) x 100cm (L)
- Easy mobility over a bed. Bassinet or examination table.
- Microprocessor Controlled
- Set Temperature range 32°C to 39.9°C.
- Temperature Alarms $\pm 0.5^{\circ}\text{C}$.
- Alarms- High/Low Temperature, Sensor Fail, Power Fail (Servo Mode)
- Programmable Features – Set Temperature (High/Low Limit), Mute, Manual Time, Apgar Time.

PHOTOTHERAPY

MODEL: IREX 12




Features

- Irradiance > 35 $\mu\text{watts/cm}^2/\text{nm}$
- Wide Effective Area
- Long Life – > 20,000 hrs.
- No UV/IR Radiation
- No Heat Generation
- Low Power Consumption
- Negligible Maintenance Cost
- Light Weight & Compact
- Height and Angle Adjustment Facility
- Mounted on Castor Wheels for Easy Mobility
- Can be used in Conjunction with any make Radiant Warmer
- Wave length : 400-550 nm
- Safety Cut Off in case of overheating (Optional)
- LCD Digital Display for Total Hours of Use and Therapeutic Hours
- Therapeutic Hours (User Resettable)
- 4 Castor Wheels



2014



OFFICE ADDRESS:

**SHOP NO.6, ALIYA MARKET, INDIRA NAGAR, NEAR AMRAPALI
CHAURAH, LUCKNOW-226016**

**Mo- 8864835791
9918555546**

PINAKI SURGICAL